

737 Pelham Boulevard, St. Paul, MN 55114 d.b.a. Memorial Blood Centers / Nebraska Community Blood Bank InnovativeBloodResources.org

Laboratory testing performed by **Memorial Blood Centers** Phone 651.332.7125 Fax 651.332.7008 **CLIA** #24D0663800

Jed Gorlin, MD, MBA, VP, Medical and Regulatory Affairs Nancy Van Buren, MD, Medical Director

REQUEST FOR REFERENCE LABORATORY TESTING

Complete information must accompany each specimen. Improperly labeled specimens will not be processed.

Date	Hospital/Lab	
Ordering Priority	City/State	
Ordering Physician	Phone () Fax ()	
TEST REQUEST		
□ Antibody Screen/Identification □ Transfusion Rea	□ NAIT Investigation (sendout)	
PATIENT INFORMATION		
Patient's Name (Last, First, Middle):		
Ethnicity: Date of Birth		
Identifying #: Sample Date, Time and Collecting Individual: (Sample noted here will be used for crossmatch, if applicable)		
Date, time, and collecting individual for each supplementary tube submitted:		
Location: ☐ Inpatient ☐ Outpatient Hgb/Hct:	Blood type:	
History of previous red cell antibody?	Antibody identified on current sample?	
Describe current transfusion problem and/or reason for submitting		
Test method(s), phase(s), and number of cells where positive reactions observed:		
Auto Control Positive? DAT: No DAT: Not performed Desitive Polyspecific Desitive IgG Desitive Complement		
CLINICAL HISTORY		
Clinical diagnosis:		
Medications:		
Has the patient been treated with immunoglobulin-based therapy in months?	ino la res (ir yes, please check one below)	
□ Anti-CD38 □ Anti-CD47 □ Anti-CD20 □ CTL □ Daratumumab □ Isatuximab □ MOR202 □ Ritu	uximab Date of Last Dose:	
Transfusion: Unknown No Yes, Date of most recent and number of RBC in past 3 months:		
Transplant: ☐ Unknown ☐ No ☐ Yes, date of transplant:	□ Autologous □ Allogeneic	
Pregnancy: ☐ No ☐ Yes, estimated date of delivery:	Number of previous pregnancies:	
Has patient received Rh immune globulin in the past 6 months? ☐ No ☐ Yes, date received:		
UNIT REQUEST		
Blood Bank unique identifer: Date and time	·	
Special Requirements: □ Irradiated □ HgbS Negative □ CMV Negative □ RH & Kell Matched □ Phenotypically Matched		
	Platelet	
☐ Crossmatched ☐ Antigen Negative	☐ Crossmatched ☐ HLA Matched	

IBR Form-0438.04 · 3/2023 Page 1 of 2

LABELING REQUIREMENTS

All samples referred for crossmatching and pretransfusion testing must meet the current Standards of the AABB regarding recipient blood samples. Sender will be notified if a sample is unacceptable; a new sample will be required.

- 1. Patient First and Last Name
- 2. Patient Identifying Number (Date of Birth not acceptable)
- 3. Date and Time Sample Collected
- 4. Phlebotomist Identity (initials)
- 5. Blood Bank Unique Identifier **crossmatch**

SPECIMEN REQUIREMENTS

Specimens collected in gel-type separation tubes are unacceptable.

Specimen may be rejected if quantity is projected to be insufficient for testing.

TEST	SAMPLE REQUIRED
ABO and Rh Typing	5-10 mL EDTA whole blood or clotted blood
Antibody Screen/Identification and Compatibility Testing	10-20 mL EDTA whole blood and 7 mL clotted blood; If patient has a positive direct antiglobulin test (DAT) include a 10-20 mL EDTA tube
Direct Antiglobulin Test	5-10 mL EDTA whole blood
Elution Study	10-20 mL EDTA whole blood
Molecular	5 mL EDTA whole blood
Hemolytic Disease of the Newborn Investigation	Mom: 10 mL clotted blood or EDTA whole blood; Baby: 2-5 mL EDTA cord blood
HLA Antigen and Antibody	10 mL EDTA whole blood and 7 mL clotted blood
Monocyte Monolayer Assay (MMA)	5 mL EDTA whole blood and two 7 mL clotted blood
Platelet Compatibility/Crossmatch	10-20 mL clotted blood or EDTA whole blood; sample must be submitted within 48 hours of collection. Samples received and frozen within 48 hours are acceptable for 7 days
Transfusion Reaction	10-20 mL clotted blood or whole blood EDTA and segments from implicated donor unit(s)

DIRECTIONS FOR SAMPLE TRANSPORT

Ship samples at ambient temperatures unless temperatures are >82F or <32F.

If ambient temperature is >82F, ship samples with coolant.

If ambient temperature is <32F, ship samples in insulated container.

STAFFING HOURS

Memorial Blood Centers' Immunohematology Reference Laboratory (IRL) is staffed from 6:00 AM Monday to 10:00 PM Friday. During these hours, contact the lab at 651.332.7125. The IRL staff is on-call from 10:00 PM to 6:00 AM Monday through Friday, weekends and holidays. For after-hours requests, call Hospital Services at 651.332.7108 and ask for the Reference Lab On-Call Technologist.

IBR Form-0438.04 · 3/2023 Page 2 of 2